

Hemostatic Dressing – QuickClot Gauze

Guideline Title: Hemostatic Dressing – QuickClot Gauze

Guideline Number: GUID-3203-PR028

Issue date: 12/26/2020	Replaces Dept. Policy:
Revision dates: 6/8/2021	Developed by: Air Care Team
Approved by: Christopher Hunter, MD Air Care Team Medical Director	Approved by:
Signature: 	Signature:
Department Numbers: 3203	

CCCXLVII. **PURPOSE:** Hemorrhage is the most common cause of preventable death in trauma. This product stops bleeding significantly faster and more effectively than standard gauze when applied with manual pressure. QuickClot products are hemostatics impregnated with kaolin, an inorganic mineral that accelerates the body's natural clotting ability.

CCCXLVIII. **DEFINITIONS:** None

CCCXLIX. **DEPARTMENT GUIDELINE:**

- A. Apply direct pressure over the wound to slow bleeding.
- B. Attempt to find the source of bleeding by wiping away any pooled blood.
- C. Open product package using the existing tear points. Keep the empty package.
- D. Begin packing wound with one end of the gauze.
- E. It is critical to pack wound deeply so the gauze makes contact with the arterial or venous breach.
- F. Do not simply cover the exterior wound with the gauze.
- G. Continue packing the wound until the entire cavity is filled.
- H. Large wound may require more than one package if QuickClot gauze.
- I. Utilize excess gauze to apply pressure directly over wound until bleeding is controlled. Use firm consistent pressure with your hand. Do not move gauze to check site. This can take 3 minutes.
- J. If the wound continues to bleed despite perceived adequate packing or if blood soaks through the gauze, remove the previous gauze placed. Wipe the area of pooled blood and find the source of bleeding. Use a new package of QuickClot and repack the wound.
- K. Once the bleeding is controlled, bandage the wound tightly to maintain pressure on the gauze and to guard against further contamination.
- L. Tuck the package of the product into the gauze to alert receiving medical personnel that the hemostatic agent has been used.

CCCL. **DOCUMENTATION:** EMS Charts

CCCLI. **REFERENCES:** <https://quikclot.com>